

Cerebral Oxygen Desaturation Events Assessed by Near-Infrared Spectroscopy During Shoulder Arthroscopy in the Beach Chair and Lateral Decubitus Positions

in JM, et al. Anesth Analg. August 2010 - Volume 111 - Issue 2 - p 496-505

Objective

The aim of this prospective cohort study was to **determine the incidence of cerebral desaturation events (CDEs) in the beach chair position (BCP)** and compare this cohort to subjects undergoing shoulder surgery in the lateral decubitus position (LDP). In addition, the relationship between CDEs and early clinical recovery was examined.

Patients

Data were collected on **124 patients undergoing elective shoulder arthroscopy** in the BCP (61 subjects) or LDP (63 subjects). Exclusion criteria included preexisting cerebrovascular disease or orthostatic hypotension, age < 18 years, ASA physical status IV or V.

Protocol

Anesthetic management was standardized in all patients. Regional cerebral tissue oxygen saturation (SctO2) was quantified using near-infrared spectroscopy (FORE-SIGHT system). Baseline heart rate, mean arterial blood pressure, arterial oxygen saturation, and **SctO2 were measured before patient positioning and then every 3 minutes for the duration of the surgical procedure. SctO2 values below a critical threshold (>20% decrease from baseline or absolute value <55% for >15 seconds) were defined as a CDE** and treated using a predetermined protocol. The number of CDEs and types of intervention used to treat low SctO2 values were recorded. The association between intraoperative CDEs and impaired postoperative recovery was also assessed.

Results

There were no differences between groups in age, weight, height, sex, preoperative hemoglobin values, preexisting medical conditions, or ASA physical status. Intraoperative hemodynamic variables did not differ between groups.

SctO2 values were lower in the BCP group throughout the intraoperative period (P < 0.0001). **The incidence of CDEs was higher in the BCP group (80.3% vs 0% LDP group)**, as was the median number of CDEs per subject (4, range 0–38 vs 0, range 0–0 LDP group, all P < 0.0001).

Table 3. Primary Outcome Variables				
	Beach chair group	Lateral group	Difference or median difference (95% CI)	P value
No. of patients	61	63	—	—
Patients with cerebral desaturation events	49 (80.3%)	0 (0%)	80.3% (68.7%–88.4%)	<0.0001
Interventions for SctO ₂ decreases	2 (0–11)	0 (0–0)	2 (2–3)	<0.0001
Interventions for MAP decreases	1 (0–6)	0 (0–9)	0 (0–1)	0.008
Episodes SctO ₂ ≤55	0 (0–4)	0 (0–0)	0 (0–0)	0.003
Episodes ≥20% decrease SctO ₂	4 (0–38)	0 (0–0)	4 (2–5)	<0.0001

CI = confidence interval; SctO₂ = regional cerebral tissue oxygen saturation; MAP = mean arterial blood pressure.
Data are number of patients (%) or median (range).

Among all study patients without interscalene blocks, **a higher incidence of nausea (50.0% vs 6.7%, P = 0.0001) and vomiting (27.3% vs 3.3%, P = 0.011) was observed in subjects with intraoperative CDEs compared with subjects without CDEs.**

Conclusions of the authors

Patients undergoing shoulder surgery in the **beach chair position may be at risk for cerebral hypoxia** because of decreases in cerebral perfusion pressure (CPP). An **association between intraoperative cerebral desaturation events and postoperative nausea and vomiting** was also observed.

Key message

Shoulder surgery in the beach chair position is associated with significant reductions in cerebral oxygenation.

For patients undergoing shoulder surgery in the beach chair position, cerebral desaturation events were observed in 80.3% of subjects.

Intraoperative cerebral desaturation events were associated with a higher incidence of nausea and vomiting in the PACU.

Pubmed Link: <https://www.ncbi.nlm.nih.gov/pubmed/20508134>

For professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult eifu.edwards.com where applicable).

Edwards devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

"ForeSight Elite tissue oximetry system is manufactured by CAS Medical, Inc., a subsidiary of Edwards Lifesciences.

Edwards, Edwards Lifesciences, the stylized E logo, ForeSight and ForeSight Elite are trademarks or service marks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2020 Edwards Lifesciences Corporation. All rights reserved. E10843/03-20/CC

Edwards Lifesciences • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com

